

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Urease from *Canavalia ensiformis* (Jack bean) – Type III, glycerol solution, 500–800 units/mL

Product Number:

U1875

CAS Number:

9002-13-5

MDL:

MFCD00070858

Storage Temperature:

2 - 8 °C

TEST

Specification

units/mL

500 - 800

One unit will liberate 1.0 micromole of ammonia from urea per minute at pH 7.0 at 25 deg C

Specification: PRD.0.ZQ5.10000067095

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.