

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calmodulin-dependent Protein Kinase II fragment 290-309 - ≥97% (HPLC)

Product Number: **C4926**
CAS Number: 115044-69-4
MDL: MFCD00076719
Formula: C₁₀₃H₁₈₅N₃₁O₂₄S
Formula Weight: 2,273.83 g/mol
Storage Temperature: -20 °C



TEST	Specification
% Peptide content	≥ 65 %
Purity (HPLC)	≥ 97 %

Specification: PRD.1.ZQ5.10000075236

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

