

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

SILu™ MAb Mouse Stable-Isotope Labeled Monoclonal Antibody – from mouse, recombinant, expressed in CHO cells

Product Number:

MSQC10

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

ug protein/vial (UV)

90 - 110

E0.1% at 280 nm = 1.4 in 0.1% formic acid

Purity by SDS-PAGE

≥ 95 %

% Incorporation

≥ 98.0 %

Specification: PRD.0.ZQ5.10000081128

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.