

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Diamond - nanopowder, <10 nm particle size (TEM), ≥95% trace metals basis

Product Number: 636444
CAS Number: 7782-40-3
MDL: MFCD00211867
Formula: C
Formula Weight: 12.01 g/mol

TEST	Specification
Appearance (Color)	Conforms to Requirements
Grey to Black	
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Trace Metal Analysis	≤ 50000.0 ppm
Purity	Meets Requirements
≥95% Based on Trace Metals Analysis	
Average Particle Size	Conforms
< 10 nm (TEM)	

Specification: PRD.0.ZQ5.10000017899

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.