

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Benzil - 98%

Product Number: B5151 CAS Number: 134-81-6

MDL: MFCD00003080 Formula: C14H10O2 Formula Weight: 210.23 g/mol

TEST Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (GC)
Solubility (Turbidity)

c = 50 mg/ml; Ethanol

Solubility (Color)

Specification: PRD.0.ZQ5.10000039931

Faint Yellow to Dark Yellow

Powder or Crystals or Granules or Chunks

Conforms to Structure

> 97.5 %

Clear to Slightly Hazy

Yellow to Very Dark Yellow and

Yellow-Green