

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Deoxyribonucleic acid, single stranded from human placenta

Product Number:

D3287

CAS Number:

9007-49-2

MDL:

MFCD00130921

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

DNA Concentration

mg/ml DNA

9 - 11 mg/ml

Hepatitis B virus

Negative

Human Source (HCV)

Negative

Human Source (HIV 1/2)

Negative

Dot Blot

Suitable

Suitable for nucleic acid hybridization.

Detects 50 pg of DNA

Specification: PRD.0.ZQ5.10000062536

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.