


Product Specification

Product Name:
4-Phenylphenol - 97%

Product Number: 134341
CAS Number: 92-69-3
Formula: C₁₂H₁₀O
Formula Weight: 170,21 g/mol



| TEST | Specification |
|---------------------|-----------------------|
| Appearance (Colour) | White to Light Brown |
| Appearance (Form) | Powder or Flakes |
| Purity (GC) | ≥ 96.5 % |
| Infrared Spectrum | Conforms to Structure |

Specification: PRE.0.ZQ5.10000021142

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

