

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fibrinopeptide A human - ≥96% (HPLC)

Product Number: F3254
CAS Number: 25422-31-5
MDL: MFCD00076437
Formula: C63H97N19O26
Formula Weight: 1,536.56 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/mL, H ₂ O	
Peptide Content	≥ 70 %
By Amino acid analysis	
Purity (HPLC)	≥ 96 %

Specification: PRD.1.ZQ5.10000038127

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.