

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Calcium thiocyanate tetrahydrate - 95%

Product Number:**520144**

CAS Number:

65114-14-9

MDL:


MFCD00661118

Formula:

C₂CaN₂S₂ · 4H₂O

Formula Weight:

228.30 g/mol

**TEST****Specification**

Appearance (Color)

White to Yellow

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Titration with EDTA

94.5 - 105.5 %

Specification: PRD.0.ZQ5.10000038623

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.