

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Diffinity RapidTip® - for PCR Purification

Product Number: **D1947**

TEST	Specification
------	---------------

dNTP Removal < 13% at 60 seconds	Conforms
Primer Removal < 13% at 60 seconds	Conforms
dsDNA Recovery > 80% at 60 seconds	Conforms

Specification: PRD.0.ZQ5.10000084465

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.