

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Anti-phospho-p53 (pSer<sup>392</sup> ) antibody produced in rabbit – affinity isolated antibody, buffered aqueous glycerol solution

Product Number: P8982

MDL: MFCD02097164

TEST	Specification
Formulation	Conforms
Solution in 10 mM Sodium HEPES, pH 7.5, 150 mM NaCl, 100 ug/mL BSA, and 50% glycerol	
Working Dilution	Conforms
by Direct Immunoblotting	
1:1,000 using p53 fusion proteins with CKII phosphorylation or an extract from hydroxyurea-treated Mv1Lu cells and chemiluminescent detection	

Specification: PRD.0.ZQ5.10000088136

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.