


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Cyclopentene-1,2-dicarboxylic anhydride - 96%

Product Number: 318353
CAS Number: 3205-94-5
MDL: MFCD00013330
Formula: C₇H₆O₃
Formula Weight: 138.12 g/mol



TEST

Specification

Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals or Flakes
Infrared spectrum	Conforms to Structure
Morpholine Titration by HCl	95.5 - 104.5 %
Purity (GC)	≥ 95.5 %

Specification: PRD.1.ZQ5.10000017595

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.