

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Perfluorinated resin solution containing Nafion™ 1100W - 5 wt. % in lower aliphatic alcohols and water, contains 15-20% water

Product Number:

274704

MDL:

MFCD00132790



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer) Approximately 15% - 25%	
Residue on Evaporation Approximately 4% - 6%	
Registered Trademark Nafion® is a registered trademark of The Chemours Company FC LLC	Confirmed

Specification: PRD.2.ZQ5.10000019877

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

