


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
6-Amino-2-naphthoic acid - 90%

Product Number: 512370
CAS Number: 116668-47-4
MDL: MFCD01861831
Formula: C₁₁H₉NO₂
Formula Weight: 187.19 g/mol



TEST

Specification

| | |
|--|--|
| Appearance (Color) | White to Beige and Faint Brown to Light Brown |
| Appearance (Form) | Powder or Crystals or Powder with Chunk(s) or Crystals with Chunk(s) |
| Infrared spectrum | Conforms to Structure |
| Purity (Titration by HClO ₄) | 89.5 - 110.5 % |
| Purity (HPLC) | ≥ 89.5 % |

Specification: PRD.0.ZQ5.10000028432

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.