

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Arsenic – 99.999% trace metals basis

Product Number: 202657
CAS Number: 7440-38-2
MDL: MFCD00085309
Formula: As
Formula Weight: 74.92 g/mol

As

TEST	Specification
------	---------------

Appearance (Color)	Light Grey to Black
Appearance (Form)	Powder or Granules or Chunk(s)
ICP Major Analysis	Confirmed
Confirms Arsenic Component	
Purity	Meets Requirements
99.999% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.1.ZQ5.10000079031

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.