

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(S)-3-(Boc-amino)-2,4-dimethyl-2-pentanol - 97%

Product Number: **719293**
CAS Number: 157035-77-3

Formula: C₁₂H₂₅NO₃
Formula Weight: 231.33 g/mol



TEST	Specification
Appearance (Color)	White to Tan
Appearance (Form) Powder or Solid	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	60.1 - 64.5 %
Nitrogen	5.8 - 6.3 %
Optical Rotation c = 0.5%, CHCl ₃	-17.2 - -11.2 deg
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000032451

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.