

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:
Coomassie protein assay reagent

Product Number: 27813

MDL: MFCD00163371

Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

Red to Very Dark Red

Appearance (Turbidity)

Clear

Appearance (Form)

Liquid

Suitability

Suitable

Suitable as protein determination reagent according to Bradford

Specification: PRD.0.ZQ5.10000037427

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.