

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Zinc nitrate hexahydrate - reagent grade, 98%

Product Number:**228737** $Zn(NO_3)_2 \cdot 6H_2O$

CAS Number:

10196-18-6

MDL:

MFCD00149889

Formula:

 $N_2O_6Zn \cdot 6H_2O$

Formula Weight:

297.49 g/mol

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Powder or Crystals or Flakes

X-Ray Diffraction

Conforms to Structure

Complexometric EDTA

21.4 - 22.5 %

% Zn

Solubility (Turbidity)

Clear to Slightly Hazy

10% in H₂O

Solubility (Color)

Colorless

Specification: PRD.1.ZQ5.1000007413

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.