

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Rabbit IgG (whole molecule) antibody produced in goat - IgG fraction of antiserum, buffered aqueous solution

Product Number:

R5506

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	8.0 - 16.0 mg/ml
Antibody Concentration(Method)	Quantitative Precipitin Assay
Antibody Concentration	3.0 - 4.5 mg/ml
ELISA - Titer	1:30,000-1:60,000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000018493

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.