

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lysine Oxidase from *Trichoderma viride* - lyophilized powder, ≥ 20 units/mg protein

Product Number:

L6150

CAS Number:

70132-14-8

TEST**Specification**

Appearance (Form)	Lyophilized Powder
Appearance (Colour)	Faint Yellow to Tan
Protein Content (Method)	Lowry
Prot. Content (% Protein)	5.00 - 20.00 %
Enzym. Assay (Spec. Activity)	≥ 20.00 U/mgP
Recommended Retest Period	-----
2 Years	

Specification: PRE.2.ZQ5.1000000928

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.