

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium - foil, thickness 0.25 mm, 99.99% trace metals basis

 Product Number:
 267481

 CAS Number:
 7440-32-6

 MDL:
 MFCD00011264

Formula: T

Formula Weight: 47.87 g/mol



TEST Specification

Appearance (Color)

Grey

Appearance (Form)

Foil

X-Ray Diffraction

Thickness

0.25 mm

Trace Metal Analysis

Purity

99.99% Based On Trace Metals Analysis

Specification: PRD.1.ZQ5.10000033469

Conforms

Conforms

Conforms to Structure

Conforms

 \leq 150.0 ppm

Meets Requirements

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.