

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Anti-Human IgG (Fc Specific)–Agarose antibody produced in goat – affinity isolated antibody, PBS suspension

Product Number: A3316

TEST	Specification
Appearance (Form)	Suspension
Appearance (Color)	White
Immobilized Protein mg/ml Resin	3.0 - 4.0
Binding Capacity mg/ml resin	2.4 - 3.4
Elution-Elution Reagent	0.2M Glycine-HCl pH 2.2
Elution-Capacity mg/ml resin	1.9 - 2.9
Recommended Retest Period 3 Years	-----

Specification: PRE.1.ZQ5.10000015976

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.