


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Oil Red O solution – 0.5% in isopropanol

Product Number: **O1391**



Formula Weight: 408.49 g/mol

TEST	Specification
Appearance	Conforms
dark red, sediment free	
Absorbance	0.20 - 0.50
diluted 0.25ml to 50 ml with Isopropyl alcohol	
Lambda Max	510 - 517 nm

Specification: PRD.0.ZQ5.10000045263

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.