

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-[4-(1H,1H,2H,2H-Perfluorodecyl)phenyl]-1,1-diphenylmethanol - 98% (HPLC)

Product Number:

672475

CAS Number:


649561-66-0

Formula:

C₂₉H₁₉F₁₇O

Formula Weight:

706.43 g/mol



TEST

Specification

Appearance (Color)

White to Beige

Appearance (Form)

Crystal to Powder

Proton NMR spectrum

Conforms to Structure

Fluorine NMR Spectra

Conforms to Structure

Purity (HPLC)

≥ 97.5 %

Supplier information

Confirmed

Product of FLUOROUS Technologies, Inc.

Specification: PRD.0.ZQ5.10000054127

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.