

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Myrcene - sum of isomers, ≥90%, natural, FG


Product Number: **W276212**

CAS Number: 123-35-3

Formula: C₁₀H₁₆

Formula Weight: 136.24 g/mol

Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Liquid
Refractive index at 20 °C	1.450 - 1.480
Infrared Spectrum	Conforms to Structure
GC (% Major 1)	≥ 88.0 %
Myrcene	
Total Purity	≥ 90.0 %
Sum of Myrcene and C ₁₀ H ₁₆ isomers	
Specific Gravity	0.785 - 0.810
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.2.ZQ5.10000100338

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.