

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fast Red Violet LB Base solution

Product Number: 912

MDL: MFCD00007780



Storage Temperature: 2 - 8 °C

TEST

Specification

Appearance (Color)	Light Yellow to Pink
Appearance (Form)	Solution
Appearance (Turbidity)	Clear
Essentially free of particles	
Wavelength	341 - 346 nm
Absorbance	1.390 - 1.580
Performance	Conforms
When staining with Sigma-Aldrich Procedure 91, sites of activity show	bright red granulation; activity weak or absent in monocytes and lymphocytes
Expiration Date	Accept

Specification: PRD.2.ZQ5.10000027840

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.