

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-LAMP1 antibody produced in rabbit – affinity isolated antibody, buffered aqueous solution

Product Number: L1418

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.0 - 1.4 mg/ml
ICC/IF Cells	NRK rat kidney cells
ICC/IF-Working Concentration	5.0 - 10.0 µg/ml
ICC/IF Cells	NIH-3T3 Mouse Fibroblasts
ICC/IF-Working Concentration	5.0 - 10.0 µg/ml
ICC/IF Cells	Hela Cells
ICC/IF-Working Concentration	5.0 - 10.0 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000015034

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.