

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Zirconium(IV) chloride – anhydrous, powder, 99.99% trace metals basis

Product Number:

647640

CAS Number:

10026-11-6

MDL:


MFCD00011306

Formula:

Cl₄Zr

Formula Weight:

233.04 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Titration by AgNO ₃ % Cl	60.2 - 61.5 %
ICP Major Analysis Confirms Zirconium Component	Confirmed
Trace Metal Analysis	≤ 150.0 ppm
Purity 99.99% Based On Trace Metals Analysis	Meets Requirements

Specification: PRD.0.ZQ5.10000036500

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.