

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Graphene dispersion - 10 mg/mL, dispersion in NMP

Product Number:

803839

Storage Temperature:

2 - 8 °C

TEST**Specification**

Appearance (Color)	Grey to Black
Appearance (Form)	Liquid or Suspension
Id/Ig Ratio (Raman) < 0.2	Confirmed
Graphene Resistance < 10 ³ ohm/sq	Confirmed
Product of Supplier Ocean NanoTech, LLC	Confirmed
Recommended Retest Period 1 Year	-----

Specification: PRD.0.ZQ5.10000072527

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.