


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lithium bis(trimethylsilyl)amide - 97%

Product Number: 324620
CAS Number: 4039-32-1
MDL: MFCD00008261
Formula: C₆H₁₈LiNSi₂
Formula Weight: 167.33 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder or Crystalline Chunks or Crystals
Titration Sec-Butanol	96.5 - 103.5 %

Specification: PRD.2.ZQ5.10000004123

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.