


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
4-Acetylphenylboronic acid

**Product Number:** 470821  
**CAS Number:** 149104-90-5  
**MDL:** MFCD01074667  
**Formula:** C<sub>8</sub>H<sub>9</sub>BO<sub>3</sub>  
**Formula Weight:** 163.97 g/mol



### TEST

### Specification

Appearance (Color)

Conforms to Requirements

White to Off-White, Beige or Light Yellow

Appearance (Form)

Conforms to Requirements

Powder and/or Chunks

Proton NMR Spectrum

Conforms to Structure

Titration with NaOH

Typically > 95%

Specification: PRD.0.ZQ5.10000037953

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.