

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Human Serum – from platelet poor human plasma, sterile-filtered, (mycoplasma tested, virus tested)

Product Number: **P2918**

MDL: MFCD00165829

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Appearance (Color)	Faint Yellow to Dark Yellow-Brown
Appearance (Form)	Liquid
pH	7.0 - 9.0
Osmolality	260 - 360 mOs/kg
Endotoxin Level	≤ 10 EU/ml
Tested For Infectious Agents	Tested
Sterility	Sterile
by DI	
Mycoplasma Test	None Detected
Cholesterol	110 - 210 mg/dl
Glucose	60 - 140 mg/dl
Iron (ug/dl)	30 - 100
Total Protein g%	4.0 - 6.5
Triglyceride	≤ 150 mg/dl
Sodium (Na)	100 - 200
Expressed in mmol/L	

Specification: PRD.1.ZQ5.10000013509

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.