

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,8-Diamino-p-menthane, mixture of cis and trans isomers - technical grade, 85%

Product Number:

D19605

CAS Number:

80-52-4

MDL:


MFCD00001488

Formula:

C₁₀H₂₂N₂

Formula Weight:

170.3 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 82.5 %

sum of cis and trans isomers

Purity (Titration by HClO₄)

82.5 - 117.5 %

Specification: PRD.1.ZQ5.10000010200

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.