

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
IRON(III) CHLORIDE SOLUTION 45 % FECL₃,

Product Number: 12322

MDL: MFCD00011005

Formula: Cl₃Fe

Formula Weight: 162.20 g/mol

| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|---|---------------|
| Titration by Na ₂ S ₂ O ₃ % FeCl ₃ | 44.6 - 45.4 % |
|---|---------------|

Specification: PRD.0.ZQ5.10000080592

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.