

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Chloromethane solution - 1.0 M in diethyl ether

Product Number:

294799



MDL:

MFCD00000872

Formula:

CH₃Cl

Formula Weight:

50.49 g/mol

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Water (by Karl Fischer)	≤ 0.005 %
Mass Balance	0.95 - 1.10 M

Specification: PRD.1.ZQ5.10000012732

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.