

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cortisone-2,2,4,6,6,9,12,12-d₈ - ≥98 atom % D, ≥98% (CP)

Product Number:

900170

Formula:


C₂₁D₈H₂₀O₅

Formula Weight:

368.43 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Form)

Powder

Purity (HPLC)

≥ 98.0 %

Proton NMR Spectrum

Conforms to Structure

D Enrichment

≥ 98.0 %

Mass Spectrometry

Conforms to Structure

Specification: PRD.0.ZQ5.10000102069

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.