

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Insulin-like Growth Factor-I Receptor human -  $\geq 95\%$  (SDS-PAGE), recombinant, expressed in NSO cells, lyophilized powderProduct Number: **I4657**

MDL: MFCD01323726

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Lyophilized from 0.2  $\mu\text{m}$  filtered solution of Phosphate Buffered Saline.

Bioactivity	Conforms
-------------	----------

Binds IGF-1 in solution.

Purity by SDS Electrophoresis	$\geq 95\%$
-------------------------------	-------------

Endotoxin Assay	Conforms
-----------------	----------

< 0.1 eu/ $\mu\text{g}$  Cytokine

Specification: PRD.0.ZQ5.10000043358

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.