

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Cytokeratin pan-FITC antibody, Mouse monoclonal – clone PCK-26, purified from hybridoma cell culture

Product Number:

F0397

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	2.0 - 5.0 mg/ml
F/P Molar Ratio	3.0 - 8.0
IHC - Tissue	Human Placenta or Human Tonsil or Human Skin
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:25-1:50
Recommended Retest Period	-----
3 Years	

Specification: PRE.1.ZQ5.10000020290

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

