


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
5 $\alpha$ -Androstan-3 $\beta$ -ol – powder

Product Number: **A2480**  
CAS Number: 1224-92-6  
MDL: MFCD00067608  
Formula: C<sub>19</sub>H<sub>32</sub>O  
Formula Weight: 276.46 g/mol



### TEST

### Specification

Appearance (Color)  
Appearance (Form)  
Solubility  
Soluble at 14 mg/mL in DMSO with warming to 60 Deg C  
Infrared spectrum  
Carbon

White to Faint Red  
Powder  
Conforms  
  
Conforms to Structure  
80.5 - 84.2 %

Specification: PRD.0.ZQ5.10000023369

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.