

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Fluorostyrene - contains 4-tert-butylcatechol as inhibitor, 98%

Product Number:**219452**

CAS Number:

350-51-6

MDL:

MFCD00000341

Formula:


C₈H₇F

Formula Weight:

122.14 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless to Faint Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

> 97.5 %

Color Test

≤ 150 APHA

Inhibitor

50 - 150 ppm

4-tert-butylcatechol

Specification: PRD.1.ZQ5.10000062907

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.