

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Propylphenol - ≥97%, FG

Product Number: W364908
CAS Number: 645-56-7
MDL: MFCD00002395
Formula: C9H12O
Formula Weight: 136.19 g/mol



TEST

Specification

Appearance (Color)	Colorless to Yellow and Colorless to Orange
(May Darken On Storage with no loss in purity)	
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 ° C	1.522 - 1.528
Purity (GC)	≥ 97.0 %
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.1.ZQ5.10000036990

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.