


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Aloin - from Curacao aloe, ~50%

Product Number: A0451
CAS Number: 1415-73-2
MDL: MFCD00151160
Formula: C₂₁H₂₂O₉
Formula Weight: 418.39 g/mol



TEST	Specification
Appearance (Color)	Yellow to Brown
Appearance (Form)	Powder
Solubility (Color)	Dark Red to Dark Purple
Solubility (Turbidity)	Clear
50 mg/ml Pyridine	
Loss on Drying	≤ 3.0 %
Barbaloin Content	40 - 60 %

Specification: PRD.0.ZQ5.10000020706

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.