

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Amino-4'-methoxyacetophenone hydrochloride - 90%, technical

Product Number:

415944

CAS Number:

3883-94-1

MDL:


MFCD00193078

Formula:

C₉H₁₁NO₂ · HCl

Formula Weight:

201.65 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Yellow or Yellow-Tan

Appearance (Form)

Conforms to Requirements

Powder or Crystalline Powder

Infrared spectrum

Conforms to Structure

Carbon

47.7 - 59.5 %

Nitrogen

6.2 - 7.7 %

Specification: PRD.0.ZQ5.10000051033

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.