

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Calcium iodide hydrate - Vetec[™] reagent grade

Product Number:

V900071

CAS Number:

71626-98-7



Formula:

CaI₂ · xH₂O

Formula Weight:

293.89 g/mol

TEST

Specification

Appearance (Colour)

White to Off-White to Tan

Conforms to Requirements

Appearance (Form)

ICP Major Analysis

Confirms Calcium Component

Compleximetric EDTA

% Ca Anhydrous Basis

Crystals or Crystals with Chunks

Confirmed

13.3 - 14.0 %

Specification: PRE.1.ZQ5.10000002586

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.