


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Aminobenzyl alcohol - 98%

Product Number: 191558
CAS Number: 623-04-1
MDL: MFCD00014782
Formula: C₇H₉NO
Formula Weight: 123.15 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White to Yellow to Brown
Appearance (Form)	Powder or Crystals or Chunk(s)
Infrared Spectrum	Conforms to Structure
Titration with HCl	97.5 - 102.5 %
GC (area %)	≥ 97.5 %

Specification: PRD.4.ZQ5.10000065671

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.