

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

S,S-DIBENZYL TRITHIOCARBONATE, 97%

Product Number: 746304 CAS Number: 26504-29-0

Formula: C15H14S3 Formula Weight: 290.47 g/mol Storage Temperature: $2 - 8 \ ^{\circ}$

TEST Specification

Appearance (Color) Appearance (Form) Infrared spectrum Carbon

Sulfur Purity (GC)

Specification: PRD.0.ZQ5.10000039396

White to Dark Yellow Powder or Crystals Conforms to Structure 59.8 - 64.2 % 32.0 - 34.3 % ≥ 96.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.