


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-(β-Chlorophenyl)ethylamine - 97%

Product Number: 340847
CAS Number: 13078-79-0
MDL: MFCD00047957
Formula: C₈H₁₀ClN
Formula Weight: 155.62 g/mol



TEST

Specification

Appearance (Color)	Colorless to Very Faint Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20 °C	1.547 - 1.551
Purity (GC)	≥ 96.5 %
Purity (Titration by HCL04)	96.5 - 103.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.