

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Glycerol - BioReagent, suitable for cell culture, suitable for insect cell culture, suitable for electrophoresis, ≥99% (GC)

**Product Number: G2025** CAS Number: 56-81-5

MDL: MFCD00004722

Formula: C3H8O3 Formula Weight: 92.09 g/mol он но Дон

TEST Specification

Appearance (Clarity) Clear
Appearance (Color) Colorless

Appearance (Form) Viscous Liquid GC (area %)  $\geq$  99 % Heavy Metals < 5 ppm

as Pb

Solubility (Turbidity) Clear

1mL/mL in H2O

Solubility (Color) Colorless

Sulfate  $\leq$  0.005 % Chlorinated Compounds < 0.005 %

Insect Cell Test Pass
Performance Pass

Suitability Suitable
Suitable for Use in Electrophoresis

Recommended Retest Period ------

3 Years

Specification: PRD.1.ZQ5.10000020951

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.