

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithocholic acid-2,2,4,4-d4 - 98 atom % D

Product Number:

589349

CAS Number:

83701-16-0

MDL:


MFCD00190473

Formula:

C24D4H36O3

Formula Weight:

380.55 g/mol



| TEST | Specification |
|-----------------------|-----------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder |
| Purity (HPLC) | ≥ 99 % |
| Proton NMR Spectrum | Conforms to Structure |
| Initial Melting Point | |
| Final Melting Point | |
| D Enrichment | ≥ 98.0 atom % |

Specification: PRD.0.ZQ5.10000102498

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

