

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Tetramethylammonium hydroxide solution - 25 wt. % in methanol

Product Number:

334901

MDL:


MFCD00008280

Formula:

C₄H₁₃NO

Formula Weight:

91.15 g/mol



TEST

Specification

Appearance (Color)

Colorless to Pale Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Titration by HCL

23.0 - 27.0 %

Water (by Karl Fischer)

≤ 10.0 %

Specification: PRD.0.ZQ5.10000013063

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.